

# Evaluation of the cleanliness of the surface by contact angle measurements

**Author** N. Bhandari

**Device** Contact Angle Meter – DMo-502

**Point** The surface evaluation is more

accessible and reliable by measuring

its water contact angle

**Keywords** Contact angle, wettability, wetting,

film thickness, surface modification, cleaning evaluation, stain evaluation



## Background

The contact angle measurement is very sensitive and affected by nanometer film thickness, and even slight changes in the surface can be detected. Therefore, it can be used as a measurement to directly investigate the affinity between liquid and solids, such as ink and the to-be-coated surface, and also an indirect measurement by determining changes in solid surfaces caused by contamination, cleaning process, surface modification treatment (corona or plasma discharge), etc. by comparing contact angles of distilled water (*Fig.* 1).




Fig. 1 – Hydrophilic by plasma treatment

### Problems to be solved

For instance, wafers and substrates must be thoroughly cleaned, or defects will occur in the following processes. Although it is possible to measure the residual oil content after cleaning by X-ray or infrared spectroscopic analysis, these devices require a reasonably skilled measurement technique, are large in scale, and very expensive.

In contrast, a simple evaluation index of surface treatment and cleanliness using a contact angle meter is relatively easy to measure and relatively inexpensive, so it is a device that we would like to propose in terms of cost performance.

#### Measurements and results

When testing, clean glass slides were left in the air (in a non-clean room) immediately after opening the package, and the contact angle was measured using the sessile drop method at each elapsed time (Fig. 2).




Fig. 2 – Contact angle  $\vartheta$ 

In a non-clean room, we used our contact angle meter model DM-500, a legacy instrument, under the following conditions:

Liquid sample:  $0.5~\mu\text{L}$  of distilled water

Solid sample: Unused glass slides from an unopened package

Room temperature: 20-27°C

Humidity: 30-56%






Fig. 3 – Relation of atmospheric exposure time and the contact angle

As a result, the contact angle value (Avg.), which was  $5.6^{\circ}$  immediately after opening the package, began to increase significantly after one day of exposure to air and became almost constant after two months (*Fig. 3*). On the other hand, the standard deviation (S.D.) increased during the period from 1 to 10 days and then decreased. This indicates that the contact angle measurement with water is helpful as a simple method to indirectly evaluate the surface condition of solids.

#### Conclusion

Our contact angle meter series is suitable for evaluating the cleanliness of the surface under actual usage conditions for best reliability and to ensure product quality.